



Implants			
Two-Piece Implants			
Item Name	Item Number	Item Image	Price
Implant 2.5 x 10	10050		\$70
Implant 2.5 x 12	10051		\$70
Implant 2.5 x 14	10052		\$70
Implant 3.1 x 10	10053		\$70
Implant 3.1 x 12	10054		\$70
Implant 3.1 x 14	10055		\$70
Implant 3.5 x 8	13508		\$137
Implant 3.5 x 10	13510		\$137
Implant 3.5 x 12	13512		\$137
Implant 3.5 x 14	13514		\$137
Implant 4.3 x 6	14306		\$137
Implant 4.3 x 7	14307		\$137
Implant 4.3 x 8	14308		\$137
Implant 4.3 x 10	14310		\$137
Implant 4.3 x 12	14312		\$137
Implant 4.3 x 14	14314		\$137
Implant 5.0 x 6	15006		\$137
Implant 5.0 x 7	15007		\$137
Implant 5.0 x 8	15008		\$137
Implant 5.0 x 10	15010		\$137
Implant 5.0 x 12	15012		\$137
Implant 5.0 x 14	15014		\$137
Implant 6.0 x 6	16006		\$137
Implant 6.0 x 7	16007		\$137
Implant 6.0 x 8	16008		\$137
One-Piece Implants			
One-Piece Implant 2.5 x 10	12500		\$49
One-Piece Implant 2.5 x 12	12502		\$49
One-Piece Implant 2.5 x 14	12504		\$49
One-Piece Implant 3.0 x 10	13000		\$49
One-Piece Implant 3.0 x 12	13002		\$49
One-Piece Implant 3.0 x 14	13004		\$49
Ball Implant 2.5 x 10	12506		\$49
Ball Implant 2.5 x 12	12507		\$49
Ball Implant 2.5 x 14	12508		\$49
Ball Implant 3.1 x 10	13006		\$49
Ball Implant 3.1 x 12	13007		\$49
Ball Implant 3.1 x 14	13008		\$49
Screw-Retained (SR) Components			
Abutment SR 2.7	10094		\$56
Abutment SR 3.5	10076		\$56
Abutment SR 20 degree	10087		\$125
Actuator	10206		\$25
Insert SR	10078		\$28
- Screw #2 sold separately	10203		\$12
Temporary Insert SR	10082		\$28
- Screw #2 sold separately	10203		\$28
Cap SR	10093		\$27
- Screw #2 sold separately	10203		\$12
Impression Coping SR	10092		\$32
- Screw #A sold separately	10205		\$19
Analog SR	10211		\$35
Placement Instrument SR 2.7S	10198		\$43
Placement Instrument SR 2.7L	10199		\$43
Placement Instrument SR 3.5S	10200		\$43
Placement Instrument SR 3.5L	10201		\$43

Non-Indexed Abutments			
Contoured Abutments			
Item Name	Item Number	Item Image	Price
Contoured Abutment 4.2 x 6.5	10010		\$56
Contoured Abutment 4.2 x 7.2	10011		\$56
Contoured Abutment 4.6 x 6.5	10013		\$56
Contoured Abutment 4.6 x 7.2	10014		\$56
Contoured Abutment 5.0 x 6.5	10016		\$56
Contoured Abutment 5.0 x 7.2	10017		\$56
Contoured Abutment 5.4 x 6.5	10019		\$56
Contoured Abutment 5.4 x 7.2	10020		\$56
Prepable Abutments			
Prepable Abutment 4.2 x 6.5	10031		\$56
Prepable Abutment 4.2 x 8.5	10032		\$56
Prepable Abutment 4.6 x 6.5	10033		\$56
Prepable Abutment 4.6 x 8.5	10034		\$56
Prepable Abutment 5.0 x 6.5	10035		\$56
Prepable Abutment 5.0 x 8.5	10036		\$56
Prepable Abutment 5.4 x 6.5	10037		\$56
Prepable Abutment 5.4 x 8.5	10038		\$56
Prepable Abutment 4.2 x 11.5	10222		\$56
Prepable Abutment 4.6 x 11.5	10223		\$56
Prepable Abutment 5.0 x 11.5	10224		\$56
Prepable Abutment 5.4 x 11.5	10225		\$56
Angled Abutments			
Angled Abutment 4.2 10°	2204		\$99
Angled Abutment 4.2 20°	2201		\$99
Angled Abutment 4.6 10°	2205		\$99
Angled Abutment 4.6 20°	2202		\$99
Angled Abutment 5.0 10°	2207		\$99
Angled Abutment 5.0 20°	2208		\$99
Healing Abutments			
Healing Abutment 4.2	64246		\$43
Healing Abutment 4.7	64746		\$43
Healing Abutment 5.1	65146		\$43
Healing Abutment 5.5	65546		\$43
Healing Abutment L/B	64244		\$43
Healing Abutment UL 4.2	10142		\$43
Healing Abutment UL 4.7	10143		\$43
Healing Abutment UL 5.1	10144		\$43
Healing Abutment UL 4.2L	10145		\$43
Healing Abutment UL 4.7L	10146		\$43
Healing Abutment UL 5.1L	10147		\$43































Indexed Abutments			
Temporary Abutments, Crown			
Item Name	Item Number	Item Image	Price
Temporary Abutment, Crown 4.2	2035		\$85
Temporary Abutment, Crown 4.7	2036		\$85
Temporary Abutments, Bridge			
Temporary Abutment, Bridge 4.2	10065		\$85
Temporary Abutment, Bridge 4.7	10066		\$85
Temporary Abutment, Bridge 5.1	10067		\$85
Temporary Abutment, Bridge 5.5	10068		\$85
CI Abutments			
CI Abutment 4.2	10101		\$105
CI Abutment 4.6	10102		\$105
CI Abutment 5.0	10103		\$105
CI Abutment 5.4	10104		\$105
CI Abutment Cement-on Sleeves			
Sleeve Straight 4.2	10125		\$13
Sleeve Straight 4.6	10126		\$13
Sleeve Straight 5.0	10127		\$13
Sleeve Straight 5.4	10128		\$13
Sleeve Angled 4.2/4.6	10129		\$13
Sleeve Angled 5.0/5.4	10130		\$13
UL Abutments			
UL Abutment 4.2 x 0.6	10132		\$105
UL Abutment 4.2 x 1.2	10215		\$105
UL Abutment 4.2 x 1.8	10216		\$105
Reamer UL 4.2 x 0.6	10148		\$57
Scan Body UL	10209		\$71
- Screw #1 sold separately	10202		\$12
Denture Abutments and Housing			
ID™ Abutment			
Immediate Denture Abutment	3003		\$49
- Screw #1 sold separately	10202		\$12
Locator Abutment			
Locator Abutment 2.5	1714		\$115
Locator Abutment 3.5	1715		\$115
Locator Abutment 4.5	1716		\$115
Locator Processing Package (2-pack)	8002		\$59
Locator Processing Package (10-pack)	8003		\$266
Ball Abutments			
Ball Abutment 2	2104		\$43
Ball Abutment 3	2105		\$43
Ball Abutment 4	2106		\$43
Ball Abutment 5	2107		\$43
Ball Abutment 10°	2108		\$76
Ball Abutment Housing			
Housing 4.2 with O-ring	10080		\$22
- used with bob-KAT implants			
Housing 4.7 with O-ring	2110		\$22
- used with Ball Abutments			
O-ring 4.2 (10 pack)	10081		\$8
O-ring 4.7 (10 pack)	2109		\$8

5 - Failure to use proper placement instrument will result in a broken instrument and unusable implant!

Instrumentation

Instrumentation Highlights

- Use 1.5 and 2.0 mm drills at 1000 rpm with copious irrigation.
- First laser mark on 1.5 and 2.0 drills is at 8 mm.
- Use Final Drills at 100 rpm without irrigation.
- First laser mark on Final Drills is at 6 mm.
- Use Reamers at 100 rpm without irrigation.
- Reamer and Healing Abutment diameters have to be the same.
- Guide Pin Laser marks start at 6 mm; marks are placed every 2 mm.
- Implant Driver M is used with a Ratchet Wrench.
- Implant Driver M is designed to withstand forces up to 100 Ncm.
- Implant Driver HP is used with a Handpiece.
- Torque higher than 25 Ncm may result in Implant Driver HP fracture.
- Placement Instrument is abutment specific. Consult Restorative Manual.
- Removal Instrument is used to remove all abutments.
- Depth Gage Laser marks start at 6 mm; marks are placed every 2 mm.
- Milling Analog should be used whenever extensive abutment preparation is required.
- Sonicate milled abutments in an alcohol bath for 3–5 minutes after each milling session.
- Implant and Milling Analogs are compatible with all abutments of KAT System.

Item Name	Item Number	Item Image	Price
Twist Drills	Specify Size		\$35
Final Drills	Specify Size		\$50
Reamer 4.2 mm	4043		\$57
Reamer 4.7 mm	4048		\$57
Reamer 5.1 mm	4052		\$57
Reamer 5.5 mm	4056		\$57
Reamer L/B	4067		\$57
Reamer UL 4.2 x 0.6	10148		\$57
Trial Implants 2.5 - 6.0 mm	Specify Size		\$42
Drill Extender	4180		\$40
Guide Pin	4300		\$16
Screw #1	10202		\$12
Screw #2	10203		\$12
Screw #A	10205		\$19
Screw #B, Short	40075		\$19
Screw #B, Long	4007L		\$19
Cover Screw	7002		\$12
Removal Screw, Short	40085		\$19
Removal Screw, Long	4008L		\$19
Bite Screw 3mm	10207		\$25
Bite Screw 6mm	10208		\$25
Open Tray Impression Coping, 7mm	10120		\$32
Open Tray Impression Coping, 12mm	10121		\$32
Screw #B sold separately			
Impression Coping IL	4015		\$28
Screw #B sold separately			
Impression Coping TT 3.5	10117		\$35
Impression Coping TT 5.5	10118		\$35
Impression Coping TT Adjustment Tool	10111		\$27
Closed Tray Impression Coping	10210		\$30
Screw #1 sold separately			
Square Driver FG	4011		\$67
Square Driver, Short	4010S		\$52
Square Driver, Long	4010L		\$52
Implant Driver M, Short	4003S		\$75
Implant Driver M, Long	4003L		\$75
Implant Driver HP, Short	4201		\$75
Implant Driver HP, Long	4203		\$75
Placement Instrument, Short	4001S		\$43
Placement Instrument, Long	4001L		\$43
Removal Instrument, Short	4002S		\$43
Removal Instrument, Long	4002L		\$43
Milling Analog	5003		\$67
Instrumentation O-ring (10)	9100		\$20
Implant Analog	5001		\$38
Tapping Instrument #1	4021		\$76
Tapping Instrument #2	8015		\$76
Abutment Remover	21059		\$72
Removal Tip #1	21058		\$49
Removal Tip #2	8022		\$49
Scan Body UL	10209		\$81



Aseptico Motors and Handpieces

Call for price.



Surgical Kit

Includes all instrumentation to place and restore implants.

9002
\$2,500



Restorative Kit

Includes:
Placement/ Removal inst. S, L
Placement/ Removal screw S, L
Square Driver, S, L

9003
\$268



Reamer Kit

Includes:
Reamers 4.2, 4.7, 5.1, 5.5, L/B
Placement/ Removal instrument S, L
Square Driver, FG
Impression TT adjustment tool

9011
\$355



Laboratory Kit

Includes:
Placement/ Removal inst. L
Placement/ Removal screw L
Square Driver, L

9007
\$150



Short Drill Kit

Includes short length drills in
1.5 - 6.0mm

9009
\$355