



**KAT Implants, LLC**

T 603.427.0084

E lab@katimplants.com

15 Rye Street, Suite 115  
Portsmouth, NH 03801

LAB ONLY

Packaged by \_\_\_\_\_

Date sent \_\_\_\_\_

Doctor \_\_\_\_\_

Patient \_\_\_\_\_

Address \_\_\_\_\_

State \_\_\_\_\_ Zip \_\_\_\_\_

Date Sent \_\_\_\_\_

Due Date \_\_\_\_\_

Allow 10 working days in-lab time from the time case is received + 2 days free return shipping.

Impression Coping TT Check List

- Clean implant post
- Place coping on implant all the way down
- Use Adjustment Tool if coping is not snug
- Use light and heavy body
- Check if coping is secure inside impression.
- If coping moves, retake impression!

Specific Instructions:

SHADE \_\_\_\_\_

TOOTH # (s) \_\_\_\_\_


Abutment size(s) \_\_\_\_\_

Qty.

- Emax
- Zirconia Full Contour
- Zirconia layered with Emax
- Custom Abutment - zirconia
- Custom Abutment - titanium
- Milled provisional crowns and bridges
- Soft tissue (no charge if it is KAT Abutment and restoration is made in our lab)
- Reduction coping (fabricated if prep adjustment is necessary due to insufficient occlusal reduction or undercuts. Will call to inform).

Torque all Indexed Abutments at **15 Ncm**

Torque all Non-Indexed Abutments at **25 Ncm**



Additional Instructions

Prices and services are subject to change without prior notice.

All rush cases must be pre-scheduled by calling 603-427-0084 before the case is shipped. Additional fee applies. Shipping supplies and ground labels are available upon request. Expediated inbound shipping is available for an additional fee, please call for specific pricing.